PREVENATM

INCISION MANAGEMENT SYSTEM

PATIENT GUIDE

You have been prescribed the Prevena™ Incision Management System to manage your incision. Please read this guide to learn more about this system and how to use it. This guide also includes important warnings for your safety. Call your doctor if you have any questions about the Prevena™ System.

If the therapy unit beeps, see the **Indicators and Alerts** section of this guide for more information.

Do not re-use any part of the Prevena™ System. The therapy unit, canister and dressing are for single use only. Re-use of any part could cause wound contamination, infection or failure to heal.

Device Description

Prevena™ 125 Therapy Unit

Caution
Light

Audio Pause
Light

Low Battery
Light


On / Off and Audio Pause Button:
Press once to mute audible alerts
for 60 minutes.


Press the On / Off button to start the 192 hour (eight day) life cycle of the therapy unit.

Prevena™ 125 Therapy Unit Carrying Case

Use the adjustable strap to wear the carrying case over your shoulder or across your chest. You can also wear the carrying case on your belt. Do not wrap the carrying case strap or dressing tubing around neck.

Battery Replacement


Push here to slide open the battery cover.

Replace with AA size batteries. Lithium batteries will last longer than regular batteries.

Prevena™ Incision Dressing with Pressure Indicator


When the pressure indicator is down, therapy is being applied to the dressing.

When the pressure indicator is up, therapy may not be at the proper pressure. See the **Indicators and Alerts** section on **page 5** for more details.

Pressure Indicator

Use Prevena™ Patch Strips™ (shown, for Peel & Place™ Dressing) and excess KCI Drape (not shown, for Customizable™ Dressing) to seal leaks.

Instructions for Use

Check the pressure indicator on the dressing regularly. See the **Prevena™ Incision Dressing with Pressure Indicator** section above for more details.

Do not turn off the therapy unit unless:

- your doctor instructs you
- you experience bleeding
- you have an allergic reaction

See the **Warnings** section of this guide (**page 4**) for more details. The dressing should only be applied or removed by or on the advice of your doctor.

Do not modify the therapy unit or dressing. Do not connect therapy unit or dressing to other devices that you might be using.

The therapy unit is a medical device that monitors pressure and is not a toy. Keep away from children, pets and pests as they can damage the dressing and therapy unit and affect performance. Keep the therapy unit free of dust and lint as they can also cause damage and affect performance.

Cleaning

If needed, clean the therapy unit and carrying case with a damp cloth and mild soap and water solution. Do not use bleach.

Bathing

- Bathing in tub is not recommended.
- Do not submerge the therapy unit or dressing.

Showering

- Quick showering is permitted if cleared with your doctor. Protect the therapy unit from direct spray of water. Avoid prolonged water contact with the therapy unit and dressing.
- When towel drying, be careful not to disrupt or damage the dressing.

Sleeping

- Place the therapy unit in a position where the tubing will not become kinked or pinched.
- Make sure the therapy unit cannot be pulled off the table or dropped on the floor.

Strenuous Activity

Stay away from strenuous activity when using Prevena™ Therapy. Ask your doctor:

- how active you can be
- when you can resume your daily routine

Contraindication and Warnings

Contraindication

The dressing contains silver. If you are sensitive or allergic to silver, call your doctor right away.

Warnings

Do not wrap carrying case strap or dressing tubing around neck.

Small Parts - Choking Hazard.

The Prevena™ System is a prescription medical device. Please read and follow all the instructions in this guide so the product can perform properly while in use.

If you need:

- Magnetic Resonance Imaging (MRI), your doctor or nurse must remove the therapy unit. The dressing can remain in place.
- Hyperbaric Oxygen Therapy (HBO), your doctor or nurse must remove the therapy unit **and** dressing.

If defibrillation is required, a doctor or nurse will decide if the dressing should be removed.

Bleeding

Some patients may be at risk of heavy bleeding, with or without using the Prevena™ System. This bleeding can be caused by the type of surgical procedure or by certain patient conditions.

If you see a large amount of blood near your incision or in the tubing or canister, leave the dressing in place. Turn off the therapy unit. Call your local emergency service right away.

Full Canister

The canister may fill up with fluid before your therapy is done. If you see a large amount of fluid in the canister, turn off the therapy unit. Call your doctor right away. The therapy unit will let you know when the canister is full.

Infection

The following symptoms may mean that your incision is infected. Call your doctor right away if you:

- have swelling at the dressing
- feel feverish
- have pus at the dressing
- feel increase in soreness at the dressing
- feel itching at the dressing
- have a bad odor at the dressing
- have redness around the dressing
- feel warmth at the dressing

Allergic Reaction


The dressing may cause an allergic reaction if you are sensitive to acrylic adhesives or to silver. The following symptoms may mean you are having an allergic reaction. Call your doctor right away if you:

- notice redness
- develop a rash
- develop itching
- notice swelling
- develop hives

If you feel short of breath, your allergic reaction may be more serious. Call your local emergency service right away and turn off the therapy unit.

Indicators and Alerts

Therapy On / Off

Green light shows therapy is on.

ON / OFF AND AUDIO PAUSE BUTTON


TURN ON - Hold two seconds

TURN OFF - Hold five seconds


AUDIO PAUSE / MUTE - Press once to mute alerts

(Audio will be muted for 60 minutes.)

Alert Lights

Dressing Pressure Indicator

Up position

Check for kinked or pinched tubing.

Check for leaks.

See Leak Alert section below.

Down position?

YES

Continue therapy. If alarming for leak, see the Leak Alert section below.

NO

See Leak Alert section below.

Leak Alert

one beep, one solid light

Check all connectors for leaks.

Make sure canister is securely locked and side tabs are flush with unit.

Press firmly around dressing edge.

Use Prevena™ Patch Strips™ or excess KCI Drape to seal leak.

Leak Alert fixed?


YES

Continue therapy.

NO

Use Audio Pause; call treating physician.

Canister Full Alert

two beeps, one solid light

Check fluid level in canister.

Full or near full?

NO

Check for kinked or pinched tubing.

YES


Turn unit off. See Therapy On / Off above.

Call treating physician.

See Warnings section on page 4.

Low Battery Alert

LOW - A slow beep and one solid light - Change batteries within six - eight hours.

CRITICAL - A rapid beep and one solid light - Change batteries right away.

Hold five seconds to TURN OFF.


Press and push to open panel.

Three AA batteries required. Lithium batteries recommended.

Press and push to close panel.

Hold two seconds to TURN ON.

System Error Alert

A repeated beep increasing in volume, and two solid lights

Hold five seconds to TURN OFF.

Hold two seconds to TURN BACK ON.

Therapy running?

YES


Continue therapy.

NO

Call treating physician.

Device Life Cycle Expired

Three lights show life cycle is up. There will be a repeated beep for 15 seconds and then the unit will shut off.

Life cycle = 192 hours (eight days)

Call treating physician.

Prevena™ 125 Therapy System

Symbols Used

Refer to Clinician Guide

Fragile

Do Not Resterilize

IP24

Ingress Protection

Date of Manufacture

Non Sterile

Single Use Only

Type BF Applied Part

Consult Instructions for Use

Keep Dry

Do not use if package is damaged or open

STERILE R

Sterile Using Radiation

Use By

EC REP

Authorized Representative in the European Community

Manufacturer

Conforms with the Medical Device Directive (93/42/EEC) and has been subject to the conformity procedures laid down in the council directive

ETL Listed
Conforms to AAMI STD ES60601-1
Certified to CSA STD C22.2 No. 60601-1.

This product is designated for separate collection at an appropriate collection point.
Do not dispose of as household waste.

Specifications

Environmental and Storage Conditions

Temperature Range: -4°F (-20°C) to 140°F (60°C)

Relative Humidity Range: 15%-95%, non-condensing

Operating Conditions

Temperature Range: 41°F (5°C) to 104°F (40°C)

Atmospheric Pressure Range: 1060 hpa (-1253 ft / -381.9 m)
(for optimum performance) to 700 hpa (9878 ft / 3010 m)

Equipment not suitable for use in presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.

The Prevena™ 125 Therapy Unit is classified as a Type BF applied part under IEC 60601-1:2006.

IP24 - Protection against solid objects greater than 12.5mm and against liquid water sprays for short periods of time.

All alerts are classified as low priority according to IEC 60601-1-8:2006.

Disposal

The therapy unit and dressing should be removed by your doctor. Do not discard these items with household trash. This could violate local laws regarding hazardous waste.

Electrical Equipment Interference

Cell phones or similar products could affect the therapy unit. Move the therapy unit away from these devices if interference is suspected.

Rx Only CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician.

Contact Information

For questions about this product, supplies, maintenance or more information about KCI products and services, please contact KCI or a KCI authorized representative, or:

In the US call 1-800-275-4524 or visit www.kci1.com

KCI USA, Inc. 12930 IH 10 West, San Antonio, Texas, 78249

Outside the US visit www.kci-medical.com

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©2015 KCI Licensing, Inc. Tutti i diritti riservati. 416358 Rev A 6 / 2015
KCI Medical Products (UK) Ltd.